



PDA10JT(-EC) Amplified MCT Detector

User Guide



















Table of Contents

Chapter 1	Warning Symbol Definitions	1
Chapter 2	Safety.....	2
Chapter 3	Description.....	3
Chapter 4	Setup.....	4
Chapter 5	Operation.....	6
	5.1. <i>Output</i>	6
	5.2. <i>Gain Adjustment</i>	6
	5.3. <i>Bandwidth Filter Adjustment</i>	7
	5.4. <i>Thermoelectric Cooler</i>	7
	5.5. <i>Light-to-Voltage Conversion</i>	7
Chapter 6	Maintenance.....	8
Chapter 7	Specifications	9
	7.1. <i>General</i>	9
	7.2. <i>Detectivity</i>	10
	7.3. <i>Filter Bandwidth</i>	11
	7.4. <i>Noise Comparison</i>	12
	7.5. <i>Wavelength Sensitivity</i>	13
Chapter 8	Mechanical Drawings.....	14
Chapter 9	Regulatory.....	15
Chapter 10	Thorlabs Worldwide Contacts.....	16

Chapter 1 Warning Symbol Definitions

Below is a list of warning symbols you may encounter in this manual or on your device.

Symbol	Description
	Direct Current
	Alternating Current
	Both Direct and Alternating Current
	Earth Ground Terminal
	Protective Conductor Terminal
	Frame or Chassis Terminal
	Equipotentiality
	On (Supply)
	Off (Supply)
	In Position of a Bi-Stable Push Control
	Out Position of a Bi-Stable Push Control
	Caution: Risk of Electric Shock
	Caution: Hot Surface
	Caution: Risk of Danger
	Warning: Laser Radiation
	Caution: Spinning Blades May Cause Harm

Chapter 2 Safety



CAUTION



Mount the unit in a position that does not block airflow to the fan and heat sink!



CAUTION



Many high speed oscilloscopes have input impedances of 50Ω . In this case, do not install a 50Ω terminator. The combined loads will equal 25Ω which could allow ~ 135 mA of output current. This will damage the output driver of the PDA10JT.



CAUTION



The PDA10JT was designed to allow maximum accessibility to the photodetector by having the front surface of the diode flush with the outside of the PDA housing. When using fiber adapters, make sure that the fiber ferrule does not crash into the detector. Failure to do so may cause damage to the diode and/or the fiber. Installing an SM1RR retaining ring (included) inside the 2" threaded coupler before installing the fiber adapter will prevent damage.

Chapter 3 Description

The PDA10JT is an amplified, TEC-cooled, switchable-gain, switchable-bandwidth, HgCdTe photoconductive detector. The detector is AC coupled and requires a chopped or pulsed input detection of light signals over a wavelength range from 2.0 to 5.4 μm . Two eight-position rotary switches allow the user to vary the gain in 6 dB steps and select low pass filter bandwidth settings from 1.25 to 160 kHz. A buffered output drives 50 Ω load impedances up to 5 volt. A constant 2.5 mA bias current is provided to the detector via a Howland current pump for improved stability and low noise operation.

The detector is mounted on a thermoelectric cooler and factory set to cool the detector to -30 °C with a thermistor providing feedback to maintain a constant temperature. This cooling provides higher detectivity (D^*), which results in a lower offset at the output and allows higher gains. It also reduces thermally generated noise. The housing is used as a heat sink and includes a fan to increase the cooling capacity. It is important to note that the cooling fan will keep the heat sink at room temperature. Without it, the heat sink will heat up causing a higher temperature drop between the detector element and the heat sink, resulting in larger TEC currents. Without the fan, the TEC current will operate at its limit (~820 mA) and the detector element will no longer be temperature stabilized. Offsets will increase and fluctuate, and output noise will increase. For best results do not block, limit airflow, or stop the cooling fan.

The detector housing has a threaded mounting aperture around the diode, which is compatible with any SM1-threaded accessory. The device ships with an SM1RR Retaining Ring allowing convenient mounting of external optics, light filters, apertures, etc. The SM1-threaded mount can be easily integrated into our cage and lens tube systems.

The PDA10JT has two 8-32 (M4) tapped holes for mounting the detector on a $\text{\O}1/2$ " optical post in one of two perpendicular directions. The detector includes a 100 - 240 V, 50 - 60 Hz power supply.

Chapter 4 Setup





Figure 1 Electrical Connections

1. Unpack the detector head.
2. (Optional) Install a Thorlabs TR-series $\text{\O}1/2$ " diameter post (not included) into one of the 8-32 (M4 on –EC version) tapped holes located on the bottom and side of the head, and mount into a PH-series Post Holder (not included).
3. Connect the power supply 4-pin plug into the power receptacle on the PDA10JT.
4. Plug the power supply into a 50 - 60 Hz, 100 - 240 VAC outlet.
5. Attach a 50 Ω coax cable (i.e., RG-58U) to the output of the PDA. When running cable lengths longer than 12", we recommend terminating the opposite end of the coax with a 50 Ω resistor (Thorlabs T4119 BCN in-line terminator) for maximum performance. Connect the remaining end to a measurement device such as an oscilloscope or high-speed DAQ card.



6. Turn on the PDA10JT using the power switch located on the top side of the detector.
7. Install any desired filters, optics, adapters, or fiber adapters to the input aperture.

**CAUTION**

The PDA10JT was designed to allow maximum accessibility to the photodetector by having the front surface of the diode flush with the outside of the PDA housing. When using fiber adapters, make sure that the fiber ferrule does not crash into the detector. Failure to do so may cause damage to the diode and/or the fiber. Installing an SM1RR retaining ring (included) inside the 2" threaded coupler before installing the fiber adapter will prevent damage.

8. Apply a light source to the detector. Adjust the gain to the desired setting. **Note:** Allow a minute for the TEC controller to stabilize the temperature. For best results, allow the unit to warm up for about 30 minutes.

Chapter 5 Operation

This device is a photoconductive detector and requires a chopped or pulsed input over a wavelength range from 2.0 to 5.4 μm . The AC-coupled amplifier circuit is designed to minimize noise.



Figure 2 Gain Control, LPF Control, and Power Selector

5.1. Output

The maximum output of the PDA10JT is 10 V for high impedance loads (i.e. $R_{LOAD} > 5 \text{ k}\Omega$) and 5 V for 50 Ω loads. Adjust the gain so that the measured signal level out of the PDA10JT is below 10 V (5 V for a 50 Ω load) to avoid saturation. If necessary, use external neutral density filters to reduce the input light level. The BNC output signal is buffered with an amplifier capable of driving 50 Ω loads. A 50 Ω series resistor is included on the output to impedance match a 50 Ω coax cable. For best performance, it is recommended operating the PDA10JT with a 50 Ω terminating load located at the end of the coax cable. While this is not necessary, it eliminates ringing and distortion due to impedance mismatches.

5.2. Gain Adjustment

The PDA10JT includes a low noise, low offset, high-gain amplifier that allows gain adjustment over a 40 dB range. The gain is adjusted by rotating the gain control knob located on the side of the unit. There are 8 gain positions incremented in 6 dB steps. To adjust the gain, follow the steps below.

1. Set the gain switch to 0 dB.
2. Turn on the light source.
3. Adjust the gain setting so the output of the detector is below the saturation level as indicated by Section 5.1 above.

5.3. Bandwidth Filter Adjustment

The PDA10JT also includes an adjustable low pass filter with settings from 1.25 kHz to 160 kHz in 8 steps. This filter allows the user to optimize the PDA10JT to operate at the lowest amount of high frequency optical and electrical noise. The filter is adjusted by rotating the filter control knob, located on the side of the unit. To adjust the filter, follow the steps below.

1. Determine the maximum bandwidth required.
2. Set the filter bandwidth switch setting just above the desired bandwidth.

Note that system signal bandwidth is 80 KHz. The 160 KHz setting is provided to eliminate high frequency noise introduced above this limit.

5.4. Thermoelectric Cooler

The thermoelectric cooler built into the detector is factory set to cool the detector to -30 °C with thermistor feedback to stabilize the temperature. The housing is used as a heat sink and a fan is used to increase the cooling capacity. It is important to note that the cooling fan will keep the heat sink at room temperature. Without it, the heat sink will heat up causing a higher temperature drop between the detector element and the heat sink, resulting in larger TEC currents. Without the fan, the TEC current will operate at its limit (~820 mA) and the detector element will no longer be temperature stabilized. Offsets will increase and fluctuate, and output noise will increase. For best results, do not block, limit airflow, or stop the cooling fan. This operation is automatic and requires no input or adjustment by the user.

5.5. Light-to-Voltage Conversion

The Spectral Responsivity, $\mathfrak{R}(\lambda)$, can be obtained from Figure 6 on page 13 to estimate the amount of output voltage to expect. The light-to-voltage conversion can be estimated by factoring the wavelength-dependent responsivity of the HgCdTe detector with the gain as shown below:

$$V_{out} (V) = Gain \frac{V}{V} * \mathfrak{R}(\lambda) \frac{V}{W} * Input Power (W)$$

For low terminating resistors, <5 k Ω or 1% error, an additional factor needs to be included in the above formula. As described above, the output includes a 50 Ω series resistor (R_S). The output load creates a voltage divider with the 50 Ω series resistor as follows:

$$Scale Factor = \frac{R_{LOAD}}{(R_{LOAD} + R_S)}$$

Where, R_{LOAD} is the terminating resistor and $R_S = 50 \Omega$. For a standard 50 Ω terminator, the gain will be scaled by $\frac{1}{2}$ as shown below:

$$Scale Factor = \frac{50 \Omega}{(50 \Omega + 50 \Omega)} = 0.5$$

$$V_{out} (V) = Gain \frac{V}{V} * \mathfrak{R}(\lambda) \frac{V}{W} * Input Power (W) * Scale Factor$$

Chapter 6 Maintenance

There are no serviceable parts in the PDA10JT detector or power supply. The housing may be cleaned by wiping with a soft damp cloth. The window of the detector should only be cleaned using isopropyl alcohol and optical grade wipes. If you suspect a problem with your PDA10JT, please contact your local Thorlabs technical support office and a member of our support team will be happy to assist you.

Chapter 7 Specifications

7.1. General

General Specifications ¹	
Optical Specifications	
Wavelength Range	2.0 - 5.4 μm
Peak Wavelength (λ_p)	4.8 μm
Peak Response	300 V/W (Typical) @ λ_p
Electrical Specifications	
Gain Adjustment Range	40 dB
Gain Steps	8
Gain Settings (dB)	0, 4, 10, 16, 22, 28, 34, 40
Filter Steps	8
Filter Settings (kHz)	1.25, 2.5, 5, 10, 20, 40, 80, 160
Output Voltage	0 - 5 V (50 Ω) 0 - 10 V (Hi-Z)
Output Impedance	50 Ω
Max Output Current	100 mA
Load Impedance	50 Ω - Hi-Z
Offset ²	20 mV (45 mV Max)
Offset Drift (40 dB)	2.7 mV / $^{\circ}\text{C}$
TEC Temperature	-30 $^{\circ}\text{C}$
Bias Current	2.5 mA (0 - 40 dB)
Physical Specifications	
Detector	MCT (HgCdTe)
Active Area	1 mm ²
Surface Depth	0.11" \pm 0.02" (2.90 \pm 0.40 mm)
Output	BNC
Detector Size	3" x 2.2" x 2.2" (76.2 mm x 55.9 mm x 55.9 mm)
Weight (Detector / Power Supply)	0.42 lbs / 0.82 lbs (191 g / 372 g)
Power Supply	27 W
Input Power	100 - 240 VAC, 50 - 60 Hz
Storage Temperature	0 to 85 $^{\circ}\text{C}$
Operating Temperature	0 to 30 $^{\circ}\text{C}$

¹ All measurements performed with a 50 Ω load unless stated otherwise. The PDA10JT has a 50 Ω series terminator resistor (i.e. in series with amplifier output). This forms a voltage divider with any load impedance (e.g. 50 Ω load divides signal in half).

² After the temperature has stabilized on all gain steps. Also note that the worst case offset is on the 40 dB gain step.

Gain (Hi-Z) ³		Low Pass Filter Bandwidth		NEP Values (@ DC - 160 kHz) ⁴	
0 dB	0.8 V/V	160 k	160 kHz	0 dB	1.43×10^{-9} W/ $\sqrt{\text{Hz}}$
4 dB	1.6 V/V	80 k	80 kHz	4 dB	7.62×10^{-10} W/ $\sqrt{\text{Hz}}$
10 dB	3.2 V/V	40 k	40 kHz	10 dB	4.05×10^{-10} W/ $\sqrt{\text{Hz}}$
16 dB	6.3 V/V	20 k	20 kHz	16 dB	2.78×10^{-10} W/ $\sqrt{\text{Hz}}$
22 dB	12.6 V/V	10 k	10 kHz	22 dB	2.15×10^{-10} W/ $\sqrt{\text{Hz}}$
28 dB	25.2 V/V	5 k	5 kHz	28 dB	2.10×10^{-10} W/ $\sqrt{\text{Hz}}$
34 dB	50.1 V/V	2.5 k	2.5 kHz	34 dB	1.97×10^{-10} W/ $\sqrt{\text{Hz}}$
40 dB	100 V/V	1.25 k	1.25 kHz	40 dB	1.84×10^{-10} W/ $\sqrt{\text{Hz}}$

7.2. Detectivity




Figure 3 PDA10JT Detectivity (D^*)

³ Gain with a 50 Ω load is one-half the hi-Z gain.

⁴ NEP values measured using a 50 Ω load and a low-pass filter setting of 160 kHz; calculated at the detector's peak responsivity wavelength.

Detectivity (D^*) is defined as is defined as:

$$D^* = \frac{\sqrt{A \cdot \Delta f}}{NEP}$$

Where A is the area of the photosensitive region of the detector, Δf is the effective noise bandwidth, and NEP is the noise equivalent power.

7.3. Filter Bandwidth




Figure 4 PDA Filter Bandwidth (50 Ω Load)

7.4. Noise Comparison

Effect of Gain and Bandwidth Settings on Noise




Figure 5 Noise Comparison at Min and Max Gain and Filter Settings

7.5. Wavelength Sensitivity

PDA10JT Wavelength Sensitivity (T = -30° C)





Figure 6 Wavelength Sensitivity Graph


Chapter 8 Mechanical Drawings



Chapter 9 Regulatory

As required by the WEEE (Waste Electrical and Electronic Equipment Directive) of the European Community and the corresponding national laws, Thorlabs offers all end users in the EC the possibility to return “end of life” units without incurring disposal charges.

- This offer is valid for Thorlabs electrical and electronic equipment:
- Sold after August 13, 2005
- Marked correspondingly with the crossed out “wheelie bin” logo (see right)
- Sold to a company or institute within the EC
- Currently owned by a company or institute within the EC
- Still complete, not disassembled and not contaminated



Wheelie Bin Logo

As the WEEE directive applies to self-contained operational electrical and electronic products, this end of life take back service does not refer to other Thorlabs products, such as:

- Pure OEM products, that means assemblies to be built into a unit by the user (e.g. OEM laser driver cards)
- Components
- Mechanics and optics
- Left over parts of units disassembled by the user (PCB's, housings etc.).

If you wish to return a Thorlabs unit for waste recovery, please contact Thorlabs or your nearest dealer for further information.

Waste Treatment is Your Own Responsibility

If you do not return an “end of life” unit to Thorlabs, you must hand it to a company specialized in waste recovery. Do not dispose of the unit in a litter bin or at a public waste disposal site.

Ecological Background

It is well known that WEEE pollutes the environment by releasing toxic products during decomposition. The aim of the European RoHS directive is to reduce the content of toxic substances in electronic products in the future.

The intent of the WEEE directive is to enforce the recycling of WEEE. A controlled recycling of end of life products will thereby avoid negative impacts on the environment.

Chapter 10 Thorlabs Worldwide Contacts

For technical support or sales inquiries, please visit us at www.thorlabs.com/contact for our most up-to-date contact information.



USA, Canada, and South America

Thorlabs, Inc.
sales@thorlabs.com
techsupport@thorlabs.com

Europe

Thorlabs GmbH
europe@thorlabs.com

France

Thorlabs SAS
sales.fr@thorlabs.com

Japan

Thorlabs Japan, Inc.
sales@thorlabs.jp

UK and Ireland

Thorlabs Ltd.
sales.uk@thorlabs.com
techsupport.uk@thorlabs.com

Scandinavia

Thorlabs Sweden AB
scandinavia@thorlabs.com

Brazil

Thorlabs Vendas de Fotônicos Ltda.
brasil@thorlabs.com

China

Thorlabs China
chinasales@thorlabs.com



THORLABS
www.thorlabs.com
